Quantec PRO

General User's Guide

TO MAKE A 'CALL'

- Press the (Call) button on a call point.
- Pull out the blue switch on the Cardiac call point.
- Press the trigger switch at the end of a tail call lead.
- Pull the cord on a ceiling pull.

This will cause the call point's confidence light to slowly flash red.

Your location, together with the 'Call' message and tone, will be shown on all relevant corridor displays.

TO 'ACCEPT' A CALL AT A CORRIDOR DISPLAY

Note: Cardiac calls can only be accepted and reset at the Cardiac call point itself.

- Note the location and the level of incoming call (i.e. 'Call', 'Help Required', 'Emergency' or 'Cardiac') from the message on the display and its accompanying tone.
- Press the (Accept) button and proceed to the relevant location.



TO PUT A CALL POINT INTO 'PRESENCE'

■ Press the (Reset) button once.

This will cause the call point's confidence light to slowly flash green. Your location, together with the 'Presence' message and tone, will be shown on all relevant corridor displays.



TO CALL FOR 'ASSISTANCE' (HELP REQUIRED) AT A CALL POINT

■ If a call has already been made or the call point is in the 'Presence' state, press the button once. If not, press the (Call) button twice.

This will cause the call point's confidence light to quickly flash red. Your location, together with the 'Help Required' message and tone, will be shown on all relevant corridor displays.



TO SIGNAL AN 'EMERGENCY' AT A CALL POINT

■ Press the (Call) and (Reset) buttons at the same time.

This will cause the call point's confidence light to quickly flash red and green. Your location, together with the 'Emergency' message and tone, will be shown on all relevant corridor displays.



TO 'RESET' A CALL POINT

Note: Reset Cardiac calls by pushing in the blue switch on the Cardiac call point.

- If the call point is in the 'Presence' state, press the (Reset) button <u>once</u>.
- If the call point is <u>not</u> in the 'Presence' state, press the (Reset) button twice. This will cancel any call from the call point and the call point's confidence light will stop flashing.



TO 'DIVERT' A CALL TO ANOTHER GROUP OF CORRIDOR DISPLAYS (IF PROGRAMMED)

When there are no calls present on the system, press the button on the corridor display. The | > Setup Divert | message appears. Press the (Accept) button.

Use the button to toggle between Divert Mode: Off and Divert Mode: On

When the display shows Divert Mode: On press the (Accept) button.



TO 'CANCEL' CALL DIVERT AT A CORRIDOR DISPLAY

Repeat the above procedure. When the display shows | Divert Mode: Off press the

Important: This guide should not be considered a definitive guide to the operation of the Quantec PRO call system, as many other permutations are possible. It is the responsibility of the installer/commissioning agent to appropriately train members of staff in the correct operation of the installed system before leaving site.

